THIS DOCUMENT IS FOR TRAINING PURPOSES ONLY. PLEASE REFER TO THE OPERATOR'S MANUAL FOR COMPLETE INSTRUCTIONS.



Dissolution Units **Mixing Training Course**

P/N 460027 Rev. D



Introduction

COURSE DESCRIPTION

This course is intended to provide dry acid mixing training to inexperienced and experienced operators of the 99 and 132 Gallon Fresenius Medical Care Dry Acid Dissolution Mixers.

TEXT AND REQUIRED SUPPLIES / EQUIPMENT

- P/N 460018 Dry Acid Dissolution Unit 132 Gallons Operators Manual
- P/N 460017 Dry Acid Dissolution Unit 99 Gallons Operators Manual
- 132 Dry Acid Dissolution Unit (Installed by FMCNA Qualified Technician.)
- 99 Dry Acid Dissolution Unit (Installed by FMCNA Qualified Technician.)
- Hydrometer
- Hydrometer Cylinder
- □ Thermometer (min. req. 25° C +/-5°C (68° to 86° F) and accuracy +/- 1°C (3.6 °F) $^{\circ}$
- Bucket/Container (approx. 3.5 gal)
- pHoenix Meter (EMD pH-indicator strips, Cat. #9590 or equivalent)
- PPE Equipment (Eye Protection, Gloves)
- 1 micron filter




Dissolution Units Overview






99 and 132 GALLONS DRY ACID DISSOLUTION UNITS CONTROL PANELS

- There are two (2) preprogrammed cycles :
 - RINSE CYCLE
 - DISSOLUTION CYCLE.
- The CONTROL PANEL will display the Dry Acid Dissolution Unit STATUS at any given time.
- The right side of the control panel displays the RINSE CYCLE
- The left side of the panel displays the DISSOLUTION CYCLE.






99 and 132 GALLONS DRY ACID DISSOLUTION UNITS EXTERNAL COMPONENTS






99 and 132 GALLONS DRY ACID DISSOLUTION UNITS INTERNAL COMPONENTS






99 and 132 GALLONS DRY ACID DISSOLUTION UNITS BASIC HYDRAULICS





Mix Process

Overview





UNIT PREPARATION Step #1


Operators Manual: Section 5.0




Unit Preparation (Mix Process Step #1)

Operators Manual Section 5.0

- Power cord is connected to 120 volts, 60
 Hz, single phase 15 amp; GFI protected circuit.
- Be certain the Dry Acid Dissolution Unit
 Drain Hose is over a floor drain and
 Transfer Hose is connected to Transfer Hose Holder.
- Purified water source is turned ON.
- Power is in the ON position
- Maximum Input Water Pressure is 60
 PSI








RINSE CYCLE Step #2


Operators Manual: Section 6.0



Operators Manual Section 6.0

- 1. Before initiating the RINSE CYCLE, the operator must ensure that:
 - Access Port Lid is IN PLACE.
 - Transfer Valve is in the CLOSED position.
 - Input water source is in the ON position.





(132 Gallon Unit)


(99 Gallon Unit)




Operators Manual Section 6.0

2. Press the Rinse side **START** button.

FILL indicator (on rinse cycle side) will illuminate and tank will automatically start to fill up with water until it reaches the 25-Gallon Sensor.









Operators Manual Section 6.0


- 3. Once the 25-Gallon Sensor is reached, the process automatically steps to **RECIRCULATE** Operation.
 - The Recirculate Operation will run for twelve (12)
 minutes on the 99 gallon Dry Acid Dissolution Unit.
 - The pump and the mixer motor will run for a ten (10) minute period for the 132 gallon Dry Acid Dissolution Unit.
 - During this time period, follow the procedure outlined to inspect the spray ball.





Operators Manual Section 6.0






99 Dry Acid Unit


Use Eye Protection 132 Dry Acid Unit



Operators Manual Section 6.0

FOR TECHNICAL SERVICE ONLY!






99 Dry Acid Unit

132 Dry Acid Unit




Operators Manual Section 6.0

4. When the **RECIRCULATE** Operation is complete the Dry Acid Dissolution Unit will switch to **DRAIN** Operation and empty rinse water from the Dry Acid Dissolution Unit to the floor drain. **Note:** This is a 10 minute timed cycle for the 132 gal. Dry Acid Dissolution Unit.



- 5. At the completion of the DRAIN Operation, the Dry Acid Dissolution Unit will refill to the 25-Gallon Sensor. The Fill Indicator Light will turn on and the RINSE Operation will start.
- 5. When the second RINSE CYCLE is finished, the Dry Acid Dissolution Unit will go to the CYCLE COMPLETE Operation. The Rinse Cycle Complete Indicator Light will turn ON and the Drain Valve will remain open, allowing any residual rinse water to go down the drain






RINSE CYCLE EXAMPLE




FILL Operation






- The tank fill to the 25 Gallon Sensor.
- The unit will automatically advance to the RECIRCULATE Operation.





Example


RECIRCULATE Operation







Example

DRAIN Operation






 Recirculate Operation will stop and DRAIN valve will open.







Example

FILL Operation





- The tank fill to the 25 Gallon Sensor.
- The unit will automatically advance to the RECIRCULATE Operation.






Example


RECIRCULATE Operation **Spray Ball** Mixer Motor DISSOLUTION RINSE CYCLE COMPLETE CYCLE COMPLETE **Mixer Shaft DRY ACID** Mixing Jets DISSOLUTION UNIT Fresenius Medical Care Propellers(2)




Example

DRAIN Operation






 Recirculate Operation will stop and DRAIN Valve will open.





Example

CYCLE COMPLETE Operation





DISSOLUTION CYCLE Step #3

Operators Manual: Section 7.0




Operators Manual Section 7.0

NOTE: A COMPLETE FULL RINSE CYCLE IS RECOMMENDED BEFORE MAKING BATCH OF CONCENTRATE.

(ONLY for the 99 Gal Unit)

NOTE: A SHORT RINSE CYCLE IS PERFORMED AS PART OF THE DISSOLUTION CYCLE: FILL OPERATION, AND DRAIN IS INITIATED AND COMPLETED WHILE THE FILL INDICATOR LIGHT IS ILLUMINATED. THIS SHORTENED RINSE CYCLE SHOULD NOT BE MISTAKEN FOR A FULL RINSE CYCLE.

The "Short Rinse Cycle" is part of the Fill Operation.






Complete Full Rinse Cycle


Operators Manual Section 7.0

1. Before initiation of Dissolution Cycle:

- Ensure the tank is empty
- The 1 micron filter is installed in the concentrate filter housing.
- Check for Propellers are attached to the Mixer Shaft (for 132 Gal. Unit).
- Ensure you have the appropriate Personal Protective Equipment donned.



A) 132 Gallon Dry Acid Dissolution Unit





B) 99 Gallon Dry Acid Dissolution Unit



Operators Manual Section 7.0

2. Ensure **Power Switch is ON** and **Water Inlet Open** then place the Dry Acid Dissolution Unit in Dissolution Cycle FILL Operation by pressing the Dissolution side **START button**.






NOTE: The 132 Dry Acid Dissolution Unit will be in the **FILL** Operation until the water reaches to the Mid-Level Sensor.

NOTE: The 99 Dry Acid Dissolution Unit will do a short Recirculation Operation with only the Fill Indicator Light illuminated then fill to the Mid-Level Sensor.




Operators Manual Section 7.0

- 3. When the Mid Water Level is reached.
 - The Water inlet valve will Close.
 - The Add Granules Light will begin to flash.
 - The Unit is in PAUSE state waiting for the operator to ADD GRANULES.



A) 132 Gallon Dry Acid Dissolution Unit



B) 99 Gallon Dry Acid Dissolution Unit



Operators Manual Section 7.0

- 4. Before adding Granules make sure to check the following:
 - Ensure that no water is leaking from underneath or any external connections to and from the Dry Acid Dissolution Unit. In addition, ensure there are no leaks at the end of the Drain Hose.
 - Remove the large access lid and ensure water has stopped at the Mid-Level Sensor.
 - Depending on the input water pressure you may have to wait several minutes to verify the water does not rise above the Mid-Level Sensor.
 - Proceed to ADD GRANULES.



DISSOLUTION CYCLE ADD GRANULES

PROCEDURE CARD



Dissolution Cycle

Add Granules


Determine how much product is required.

FMCNA Dissolution Unit	No. of Cases Needed			
99 Dry Acid Dissolution Unit	6			
132 Dry Acid Dissolution Unit	8			

 Check case labels to ensure all cases have the same Catalog Number and Potassium Number.

Recommendation:

- 1. Group the cases to be used.
- 2. Separate/isolate the group of cases from other cases that are present to avoid mixing Catalog No./ Products.





Dissolution Cycle

Add Granules









.....

Dissolution Cycle

Add Granules

Fill in the <u>Production Record Form</u>, attached to the Operators Manual.

Depending on product, GranuFlo® or Citrasate® Dry Form.

Record The Following Information Depending on the Product.

					`					
	Dialysis Unit Name & Location #:				Dry Acid Dissolution Unit Serial #:					
	DRY ACID PRODUCT CASE INFO	RMATION								
	OPERATOR (print name):		DATE:						Batch #:	
	CASE 1: Dry Acid Catalog # (labe	el on box)	BOX Lo	t#				4 TH bag used √ box Yes □	Potassium # (1K, 2K, 3K)	
	CASE 2: Dry Acid Catalog #		BOX Lot #					4 TH bag used ✓ box Yes □	Potassium #	
	CASE 3: Dry Acid Catalog # CASE 4: Dry Acid Catalog #			BOX Lot#				4 TH bag used √ box Yes □	Potassium #	
				BOX Lot#				4 TH bag used ✓ box Yes □	Potassium #	
	CASE 5: Dry Acid Catalog #		BOX Lo	t#				4 TH bag used ✓ box Yes □	Potassium #	
	CASE 6: Dry Acid Catalog #		BOX Lo	BOX Lot #				4 TH bag used ✓ box Yes □	Potassium #	
	1. After Final Fill Level has been	reached, CLOSE water suppl	y valve to	the Unit. Once	e this is done check th	ne box.				
	SPECIFIC GRAVITY									
	Measured Temp	Print Catalog # -Specific Gravity	/ Value for	the Measured Te	emp listed in Appendix B:		Measured Specific Gravity Value:		Check one	
	TEMP:	LOW:		HIGH:					□-Pass	
									☐-Fail (void section)	
	OPERATOR SIGNATURE:				VERIFIER SIGNATURE:					



Add Granules

PPE

(Personal Protective Equipment is required/ recommended)



Recommended



For more information refer to the MSDS form of the Dry Acid Product.



Operators Manual Section 7.0

- A) Open a case of Citrasate® DRY / GranuFlo®.
- B) Cut off the top of 1 bag just below the bag seal.
- C) Remove the Dry Acid Dissolution Unit Small Access Lid and slowly add granules as seen in Figure.
- D) Before proceeding with the next case, repeat the steps B and C until all the bags of the case have been added.

Remember:

- Use 6 cases for the 99 Gallon Unit
- Use 8 cases for the 132 Gallon Unit All cases must have same FMCNA Catalog Number.




Figure. Add Granules



Operators Manual Section 7.0

5. On the 132 Gal Dry Acid Dissolution Unit, After the Dry Acid Concentrate has been added, dry off the Upper Level Sensor of any splashing that may have occurred during the add granules process (Figure 10). Then reinstall the Small Access Lid and press the Dissolution START button.

On the 99 Gal Dry Acid Dissolution Unit, After adding all bags of Dry Acid, reinstall Small Access Lid, then press the Dissolution START button.

Both Dry Acid Dissolution Units will proceed to the **MIX** Operation.






Figure 10. Clean Sensor






Operators Manual Section 7.0

6. For the 132 Gallon Dry Acid Dissolution Unit, during the MIX Operation, the solution is mixed for a period of forty-five (45) minutes allowing the granules to completely dissolve. The Dry Acid Dissolution Unit will then automatically step to the DEAERATION Operation.




For the **99 Gallon Dry Acid Dissolution Unit**, during the **MIX** Operation, the solution is mixed for a period of thirty-five **(35) minutes** allowing the granules to completely dissolve. The Dry Acid Dissolution Unit will then automatically step to the DEAERATION Operation.






Operators Manual Section 7.0

7. For the 132 Gallon Dry Acid Dissolution Unit, the DEAERATION Operation runs for five (5) minutes during which the entrapped air is allowed to separate out of the solution. Upon completion, the Dry Acid Dissolution Unit will automatically step to FINAL FILL Operation and the Final Fill Indicator Light will illuminate.




For the 99 Gallon Dry Acid Dissolution Unit, the DEAERATION Operation runs for two (2) minutes during which the entrapped air is allowed to separate out of the solution. Upon completion, the Dry Acid Dissolution Unit will automatically step to FINAL FILL Operation and the Final Fill Indicator Light will illuminate.





Operators Manual Section 7.0


8. In the **FINAL FILL** Operation, the supply water valve will open and Dry Acid Dissolution Unit to the final water level sensor or upper level sensor.






Operators Manual Section 7.0

9. When the final water level is reached, the unit will automatically step to HOMOGENIZE Operation. Remove the Large Access Lid and ensure the solution level has reached the Final Fill Sensor or Upper Level Sensor. Place Large Access Lid onto the Dry Acid Dissolution Unit. Then, CLOSE the water supply valve to the Dry Acid Dissolution Unit.



10. During **HOMOGENIZE** Operation, the mixer motor will stir the solution for ten (10) minutes. When the HOMOGENIZE Operation is complete, the Transfer Indicator Light will flash. Remove Large Access Lid and look into the tank to make sure the granules have dissolved and the solution is colorless. Once you have verified the granules are dissolved and colorless, the solution is ready for Specific Gravity Testing.







Dissolution Cycle (Mix Process – Step#3) **Example 132 - Gal**



Example 132 - Gallon

FILL Operation






- The tank will fill to the Mid Level Sensor.
- The unit will automatically advance to ADD GRANULES.



Example 132 - Gallon

ADD GRANULES Operation

 Remove the Small Dissolution Tank Lid and slowly add granules.






DISSOLUTION


FILL
RECIRCULATE
DRAIN
CYCLE COMPLETE

DRY ACID
DISSOLUTION UNIT
132 Gallons

Example 132 - Gallon








- Press the START button.
- The mix motor will run for 45 minutes.


NOTE: In the 99 Dry Acid Dissolution Unit the Mix Operation will run for 35 minutes.



Example 132 - Gallon

DEAERATION Operation






- DEAERATION Operation
- During this (5) minute Operation the Mixer Motor stops to allow air bubbles to separate out of the solution.


NOTE: In the 99 Dry Acid Dissolution Unit the DEAERATION Operation will run for 2 minutes.



Example 132 - Gallon

FINAL FILL Operation







- FINAL FILL Operation
- During this Operation, the Dissolution Tank is filled with treated water up to the Final Fill Sensor. Then, the unit automatically advance to the HOMOGENIZE Operation.



Example 132 - Gallon

HOMOGENIZE Operation







- HOMOGENIZE Operation
- During this ten (10) minute
 Operation, the Mixer Motor stirs
 the solution.
- At the end of this Operation, the unit will automatically advance to TRANSFER Operation.



Example 132 - Gallon

TRANSFER Operation





- TRANSFER light flashing on the front panel.
- The Solution in the tank is ready for Specific Gravity Testing.

SPECIFIC GRAVITY TEST

Operators Manual: Section 7.1



Operators Manual Section 7.1

Once the Transfer Indicator light flashes, the concentrate is ready to be tested. This test is done using a hydrometer which measures the density of the acid in the concentrate. The resultant value is called the "specific gravity". The specific gravity is measured to verify that the concentrate has been properly mixed.

Requirements for the test:

- Hydrometer
- Hydrometer Cylinder
- Thermometer*
- Bucket/Container (approx. 3.5 gal)
- •pHoenix Meter (EMD pH-indicator strips, Cat. #9590 or equivalent)
- *Minimum Requirements: 25°C \pm 5 °C (68 ° to 86 °F) and accuracy \pm 1 °C (3.6 °F)




Operators Manual Section 7.1

- 1. Prepare the Hydrometer and the Hydrometer Cylinder
 - Check Hydrometer for cracks.
 - Do not use a cracked Hydrometer.
 - Always handle the Hydrometer with care. Always hold it vertically by the top, as finger marks lower down can affect the accuracy of the instrument.
 - The Hydrometer should never be held by the stem horizontally.
 - Rinse the hydrometer and the hydrometer cylinder separately with purified water before checking solution for correct density. The purified water source must meet ANSI/AAMI or ISO standards for dialysis currently ANSI/AAMI RD62, or ISO 13959.
 - Once Hydrometer is clean, place the hydrometer onto a clean area.
 - The Hydrometer is fragile and can break easily.
 - Store in a Safe Area.



Operators Manual Section 7.1


2. Remove Transfer Line from the Transfer Hose Holder and connect the Transfer Line to the Transfer Nozzle.






Operators Manual Section 7.1

3. Hold the Transfer Nozzle and make sure the Ball Valve on the Transfer Nozzle is CLOSED.




- 4. Slightly OPEN the Ball Valve next to the Filter Housing.
- 5. Press Dissolution START button and open slightly the Transfer Nozzle Ball Valve flush approximately 3.5 gallons of solution using a container(s) that will allow for a 3.5 gallon measurement.
 - This removes any solution left in the Transfer Line from a prior batch.
 - Once the 3.5 gallons of solution have been flushed out of the transfer hose, CLOSE the Transfer Nozzle ball valve.






Operators Manual Section 7.1

- 6. Fill the Hydrometer Cylinder:
 - Insert the Transfer Nozzle into the Hydrometer Cylinder.
 - Slowly OPEN the Transfer Nozzle Valve until solution starts to slowly accumulate into the Hydrometer Cylinder. Allowing solution to fill down the side of the beaker minimizes the creation of bubbles within the solution (see figure item A).
 - Fill the Hydrometer Cylinder approximately 2/3rd full.
 - Close Transfer Nozzle and the Ball Valve next to the filter housing, then press PAUSE on the Display Panel.
 - Place Transfer Nozzle onto a clean surface.
- 7. Make sure the solution does not have excessive amount of bubbles. To release excess amount of bubbles from the solution gently tap the Hydrometer Cylinder or gently swirl (See Adjacent Picture).






Operators Manual Section 7.1

- 8. Place Hydrometer Cylinder on a level table
 - Measure the temperature of the solution in the Hydrometer and record it on the batch production record form.
 - For thermometer requirements refer to appendix A or B.
 - Go to Appendix A (GranuFlo®) or Appendix B (Citrasate® Dry).

Dialysis Unit Name & Location #:				D	ry Acid Dissolution	Uni	t Serial #:			
DRY ACID PRODUCT CASE INFO	RMATION									•
OPERATOR (print name):		DATE:			TIME:				Batch #:	
CASE 1: Dry Acid Catalog # (labe	el on box)	BOX Lo	t#					4 TH bag used √ box Yes □	Potassium # (1K	, 2K, 3K)
CASE 2: Dry Acid Catalog #			t#					4 TH bag used √ box Yes □	Potassium #	
CASE 3: Dry Acid Catalog #			t#	_			4 TH bag used √ box Yes □	Potassium #		
CASE 4: Dry Acid Catalog #			t#					4 TH bag used √ box Yes □	Potassium #	
CASE 5: Dry Acid Catalog #		BOX Lot#						4 TH bag used √ box Yes □	Potassium #	
CASE 6: Dry Acid Catalog #			t#					4 TH bag used √ box Yes □	Potassium #	
1. After Final Fill Level has been	reached, CLOSE water supply	y valve to	the Unit. Onc	e th	is is done check th	ne bo	X.			
			SPE	CIF	IC GRAVITY					
Measured Temp	Print Catalog # -Specific Gravity Value for the Measured			emp listed in Appendix B: Measured			Measured	Specific Gravity Value:	Check one	
TEMP:	LOW:		HIGH:						□-Pass □-Fail (void se	ection)
OPERATOR SIGNATURE:	•				VERIFIER SIGNAT	URE	:		•	






Operators Manual Section 7.1

9. On Appendix A or B

- Locate the Catalogue number of the dry acid product being used (1)
- Locate the value of the "measured temperature" of the solution (2)
- Moving across and then down, identify the 'low' and 'high' specific gravity numbers.
- These 'low' and 'high' numbers can be recorded in the appropriate box in the Batch Production Record Form located at the last pages on the manual (3).




Dialysis Unit Name & Location	Dry Acid Dissolu	Dry Acid Dissolution Unit Serial #:						
DRY ACID PRODUCT CASE INFO	RMATION							
OPERATOR (print name):	OPERATOR (print name):			TIME:			Batch #:	
CASE 1: Dry Acid Catalog # (lab	CASE 1: Dry Acid Catalog # (label on box)			•		4 TH bag used ✓ box Yes □	Potassium # (1K, 2K, 3K)	
CASE 2: Dry Acid Catalog #	BOX Lot #			4 TH bag used ✓ box Yes □	Potassium #			
CASE 3: Dry Acid Catalog #	BOX Lot #			4 TH bag used ✓ box Yes □		Potassium #		
CASE 4: Dry Acid Catalog #	BOX Lot #				4 TH bag used ✓ box Yes □	Potassium #		
CASE 5: Dry Acid Catalog #	CASE 5: Dry Acid Catalog #					4 TH bag used ✓ box Yes □	Potassium #	
CASE 6: Dry Acid Catalog #	CASE 6: Dry Acid Catalog #					4 TH bag used ✓ box Yes □	Potassium #	
1. After Final Fill Level has been	reached, CLOSE water suppl	y valve to the Unit. O	nce this is done che	ck the b	ox.			
			PECIFIC GRAVITY					
Measured Temp	Measured Temp Print Catalog #-Specific Gravity Value for the Measured			Temp listed in Appendix E Measured Sp			Check one	
TEMP:	LOW:	HIGH:		ı			□-Pass □-Fail (void section)	
OPERATOR SIGNATURE:	•		VERITIER SIG	TURE			•	

Minimal Thermometer Specifications: Temperature Range 25° C +/-5° C (68° to 86° F) and accuracy +/- 1° C (3.6° F)



Operators Manual Section 7.1

- 10. Obtain the clean hydrometer and slowly lower the hydrometer into the solution until the hydrometer begins to float freely.
- 11. Hold the top of the hydrometer between your thumb and finger and spin the hydrometer in the Hydrometer Cylinder.
 - The slow spinning action of the hydrometer will cause the hydrometer to stay away form the sides of the Hydrometer Cylinder.
 - In addition, this will help keep bubbles from forming on the hydrometer.






Operators Manual Section 7.1

12. Hydrometer will move up and down

If during the up and down movement of the hydrometer the meniscus is crinkled or dragged out of shape by the motion of the hydrometer, then this indicates that either the hydrometer or the surface of the solution is not clean. Clean the hydrometer and Hydrometer Cylinder again. Then restart the Specific Gravity from step 1.


 If the meniscus remains unchanged as the hydrometer rise and falls, then the hydrometer and liquid surface are clean and a reading can be taken.






Operators Manual Section 7.1

- 13. Allow the hydrometer to settle within the solution. Once the hydrometer is stable, place your head at eye level to the meniscus of the solution as seen in the adjacent Figure.
 - The point where the bottom of the meniscus crosses the hydrometer is the correct measuring point.
 - Do not take a reading if the hydrometer is touching the side of the hydrometer cylinder.






How to read the Hydrometer (SCALE)





How to read the Hydrometer

(SCALE)



MEDICAL CARE

Operators Manual Section 7.1

14. Document the Specific Gravity in the Batch Production Record Form in the "Measured Specific Gravity Value" box.

Dialysis Unit Name & Location		Dry Acid Dissolution Unit Serial #:					
DRY ACID PRODUCT CASE INFO	RMATION		<u>'</u>				
OPERATOR (print name):		DATE:			TIME:		Batch #:
CASE 1: Dry Acid Catalog # (label on box)		BOX Lo	ot#			4 TH bag used ✓ box Yes □	Potassium # (1K, 2K, 3K)
CASE 2: Dry Acid Catalog #			ot#			4 TH bag used √ box Yes □	Potassium #
CASE 3: Dry Acid Catalog #			ot#			4 TH bag used √ box Yes □	Potassium #
CASE 4: Dry Acid Catalog #			ot#			4 TH bag used √ box Yes □	Potassium #
CASE 5: Dry Acid Catalog #			ot#			4 TH bag used √ box Yes □	Potassium #
CASE 6: Dry Acid Catalog #		BOX Lo	ot#			4 TH bag used √ box Yes □	Potassium #
1. After Final Fill Level has been	reached, CLOSE water supp	ly valve to	the Unit. Once	this is done check the	e box.		
			SPEC	FIC GRAVITY	-		
Measured Temp	Measured Temp Print Catalog #-Specific Gravity Value for the Measured Ten				Measured	Specific Gravity Value:	Check one
TEMP:	LOW:		HIGH:				□-Pass □-Fail (void section)
OPERATOR SIGNATURE:				VERIFIER SIGNATU	JRE:		



Operators Manual Section 7.1

15. If the results from the Specific Gravity Test are **acceptable** (Measured Value within Low and High Values), check the "pass" box on the Dry Acid Batch Production Record Form. The solution is ready for the TRANSFER Operation.

Dialysis Unit Name & Location #:					Dry Acid Dissolution Unit Serial #:					
DRY ACID PRODUCT CASE INFO	RMATION									
OPERATOR (print name):		DATE:						tch #:		
CASE 1: Dry Acid Catalog # (label on box)		BOX Lot #				4 TH bag used ✓ box Yes □	Po	tassium # (1K, 2K, 3K)		
CASE 2: Dry Acid Catalog #			ot#				4 TH bag used ✓ box Yes □	Po	tassium #	
CASE 3: Dry Acid Catalog #			BOX Lot #				4 TH bag used ✓ box Yes □		tassium #	
CASE 4: Dry Acid Catalog #			ot#				4 TH bag used √ box Yes □	Po	tassium #	
CASE 5: Dry Acid Catalog #			ot#				4 TH bag used ✓ box Yes □	Po	tassium #	
CASE 6: Dry Acid Catalog #		BOX Lot #					4 TH bag used ✓ box Yes □	Po	tassium #	
1. After Final Fill Level has been	reached, CLOSE water suppl	y valve to	the Unit. Once	this	is done check the	box.				
			SPEC	CIFIC	GRAVITY					
Measured Temp Print Catalog #-Specific Gravity Value for the Measured Te				emp listed in Appendix B: Measured S			Specific Gravity Value:		Check one	
TEMP:	LOW:	HIGH:							□-Pass □-Fail (void section)	
OPERATOR SIGNATURE:				١	VERIFIER SIGNATU	IRE:		•		



Operators Manual Section 7.1

16. Pour residual solution from hydrometer cylinder into the solution present in the **residual solution bucket** .

Rinse the hydrometer cylinder and the hydrometer before storing equipment. To dispose the solution in the Residual Solution Bucket, See Section 13.1: Residual Solution Bucket Disposal of the Operators Manual.



General Instructions: Read carefully and follow the instructions for this exercise provided on the bottom of this page.

The Operator "X" from the Dialysis Unit "FMCNA #13" located in "Dallas, TX." is preparing to make his first batch of Citrasate® Dry Product Code OFD1231-DA using a 99 Gallon Dry Acid Dissolution Unit (mixer) with the serial number DA99-123456. The user has several cases of the product with the label shown on the next slide.

Each case of Citrasate® Dry contains 4 bags.


It is November 12, 2011 and the time is 9:25 am and the user starts with the mixing process. After mixing, the operator proceeds to test for Specific Gravity.

The results are shown below:

Temperature of Sample: 27° C, Specific Gravity: 1.198

Instructions: Using the information provided above and the label shown in the next slide, fill in the batch production record attached to the manual.







PRACTICE EXERCISE RECORD FORM

Dialysis Unit Name & Location	Гх.	Dry Acid Dissolution Unit Serial #: DA99-123456.							
DRY ACID PRODUCT CASE IN F	ORMATION		'						
OPERATOR (print name): Operator X		DATE:	11/12/201	11 TIME: 9:2			5 a.m.	Batch#: 1	
CASE 1: Dry Acid Catalog# (lal		BOX Lo		10 15 05 000			4 [™] bagused 🗸 box	Potassium # (:	1k, 2k, 3k)
	0FD1231-DA		Lot#	12ATGF009			Yes 🍑		1K
CASE 2: Dry Acid Catalog#	0ED1001 D 1	BOX Lo	t# T 04#	12 ATCE000			4 [™] bagused ✓ box	Potassium#	
	0FD1231-DA		LOI#	# 12ATGF009			Yes 🗹	1K	
CASE 3: Dry Acid Catalog#			Lot# 12ATGF009			4 [™] bagused ✓ box		Potassium#	
	0FD1231-DA		L0t#	12A1GF009			Yes 🗹		1K
CASE 4: Dry Acid Catalog#	0FD1231-DA	BOX Lo	t# I 0##	12ATGF009			4 [™] bagused ✓ box	Potassium#	
	01/D1231-DA		L0t#	12A1 OF 009			Yes 🗹		1K
CASE 5: Dry Acid Catalog#	0FD1231-DA	BOX Lo	t# T o##	12ATGF009			4 [™] bagused ✓ box	Potassium#	
	0FD1231-DA		L0t#	12A1GF009			Yes 🗹		1K
CASE 6: Dry Acid Catalog#	0FD1231-DA	BOX Lo	BOX Lot # Lot# 12ATGF009			4™ bagused 🗸 box		Potassium#	
	01D1231-DA		LOI#	12A1GF009			Yes 🗹	1K	
1. After Final Fill Level has bee	n reached, turned OFF water	valve to t	he Unit. Once th	is is done check the	box				[
			SPEC	IFIC GRAVITY					
Measured Temp	Print Catalog # - Specific Grav	rint Catalog #-Specific Gravity Value for the Measured To			mp listed in Appendix A: Measured :			Check one	
TEMP: 27° C LOW: 1.18		HIGH:		1.400			1 100	▼ -Pass	
21 0	1.187			1.199			1.198	□-Fail (void	section)
operator signature: Verifier Y Verifier Y									



APPENDIX B: CITRASATE® DRY SPECIFIC GRAVITIES TABLE

Specific Gravity Ranges

1:44 PROPORTIONING

		16.5°C to (61.7°F to		21.5°C to (70.7°F to		26.5°C to 31.4°C (79.7°F to 88.5°F)		
Catalog Number		Low	High	Low	High	Low	High	
	0FD1231-DA	1.192	1.204	1.189	1.201	1.187	1.199	
1	0FD1251-DA	1.192	1.204	1.190	1.202	1.188	1.200	
	0FD2231-DA	1.194	1.206	1.191	1.203	1.189	1.201	
2K	0FD2251-DA	1.194	1.206	1.192	1.204	1.189	1.201	
7	0FD2301-DA	1.195	1.207	1.193	1.205	1.190	1.202	
	0FD3231-DA	1.195	1.207	1.193	1.205	1.191	1.203	
3	0FD3251-DA	1.196	1.208	1.194	1.206	1.191	1.203	
m	0FD3301-DA	1.197	1.209	1.194	1.206	1.192	1.204	

Minimum Thermometer Specifications: Temperature Range 25° C +/-5°C (68° to 86° F) and accuracy +/
1°C (3.6 °F)

What thermometer to use?




TRANSFER Step # 4

Operators Manual: Section 7.2, 7.3



Transfer (Mix Process – Step #4)

Operators Manual Section 7.2 and 7.3






Transfer (Mix Process – Step #4)

Operators Manual Section 7.2 and 7.3

TRANSFER TO STORAGE TANK (SECTION 7.2, STEPS 1 TO 3)

- After you have obtained a valid specific gravity reading, remove the Transfer Line from the Transfer Nozzle.
- 1st Connect Transfer Line to the storage tank, 2nd open ball valve at the side of the Filter Housing. 3rd Press the Dissolution START button to transfer concentrate into the storage tank.
- Once Transfer of solution is complete, place Transfer line on to the Dry Acid Dissolution Unit Transfer Line Holder. Do not leave Transfer line attached to the Storage Tank.






Transfer (Mix Process – Step #4)

Operators Manual Section 7.2 and 7.3

TRANSFER TO INDIVIDUAL CONTAINER (SECTION 7.3)

- Containers must be properly labeled.
- Connect the Transfer Line to the Transfer Nozzle.
- Place the Transfer Nozzle into the opening of an individual container.
- Slightly OPEN the Ball Valve on the top of the Filter Housing.
- With the Transfer light flashing, Press Dissolution START button. The Transfer Pump will start.
- Slowly OPEN Ball Valve on the Transfer Nozzle until the desired rate of flow through the nozzle is achieved.






Transfer (Mix Process – Step #4)


Operators Manual Section 7.2 and 7.3

TRANSFER TO INDIVIDUAL CONTAINER (SECTION 7.3)

 Once the manual TRANSFER Operation has been completed and the individual containers are filled, CLOSE Transfer Nozzle Valve. Press the PAUSE button.



 When the Dry Acid Dissolution Unit is empty, the Dry Acid Dissolution Unit will step to the CYCLE COMPLETE Operation.






Maintenance Program

Operators Manual: Section 8.0



Maintenance Overview





Cleaning and Visual Inspection

Operators Manual: Section 8.1, 8.2



Maintenance

Cleaning and Visual Inspection

CLEANING

Clean the exterior surface of the Dry Acid Dissolution Unit thoroughly after each batch of concentrate is mixed. If necessary, a mild detergent solution may be used to clean the exterior surface. Care should be taken not to contaminate the system interior. All spills should be wiped off immediately. Spillage at the control panel should be avoided in order to minimize the possibility of electrical malfunction.

VISUAL INSPECTION

Visually inspect the Dry Acid Dissolution Unit prior to each use.

The operator should look for any defects which may inhibit the safe or proper Operation of the Unit. Items such as:

- Damaged hydraulic hoses or fittings.
- Damaged electrical cables or connections.
- Loose, missing, or damaged hardware.
- Previous process contamination should be corrected prior to the use of the Dry Acid Dissolution Unit.



Maintenance

Routine Maintenance Schedule

132 Dry Acid Dissolution Unit

1)2 bi y ficia bissolation onit									
PROCEDURE	PER BATCH	MONTH	AS NEEDED	REF. SECTION					
RINSE CYCLE	x			6					
VISUAL INSPECTION	x			8.1					
CLEANING SURFACES	x			8.2					
DISINFECTANT			x	8.3					
FILTER			X*	8.4					
CORROSION		X**		N/A					
SPRAY BALL			x	SECTION 6 STEP3					

^{*} It is recommended to change the Filter after mixing 6 batches or when the 132 Gallon Dry Acid Dissolution Unit Requires Disinfection. If a tank becomes contaminated, the tank will need to be disinfected before a new filter is installed.

99 Dry Acid Dissolution Unit

99 Dry Acid Dissolution Offic									
PROCEDURE	PER BATCH	MONTH	AS NEEDED	REF. SECTION					
RINSE CYCLE	x			6					
VISUAL INSPECTION	x			8.1					
CLEANING SURFACES	x			8.2					
DISINFECTANT			x	8.3					
FILTER			X*	8.6					
STAND PIPE FILTER			x	8.4.2					
CORROSION		X**		N/A					
SPRAYBALL			x	SECTION 6 STEP3					

^{*} It is recommended to change the Filter after mixing 6 batches or when the 99 Gallon Dry Acid Dissolution Unit Requires Disinfection. If a tank becomes contaminated, the tank will need to be disinfected before a new filter is installed.



^{**} It is recommended that you look for corrosion or salt deposits at the Final Fill Sensor, Propellers and Shaft within the Unit's Tank. Also, look for any corrosion around the connectors at every valve. Any excessive corroded part on the unit should be clean and replaced if needed

^{**} It is recommended that you look for corrosion or salt deposits at the Final Fill Sensor and within the Unit's Tank. Also, look for any corrosion around the connectors at every valve. Any excessive corroded part on the unit should be clean and replaced if needed.


Sodium Hypochlorite "Bleach" DISINFECTION

Operators Manual Section 8.3



Operators Manual Section 8.3

SODIUM HYPOCHLORITE (BLEACH) DISINFECTION






WARNING! ENSURE THE TRANSFER LINE IS NOT CONNECTED TO A CONCENTRATE STORAGE CONTAINER/TANK.

NOTE: USE ONLY SODIUM HYPOCHLORITE (5% TO 10%) TO DISINFECT THE GranuFlo DISSOLUTION UNIT. MAKE SURE THE BLEACH DOES NOT CONTAIN A CLEANER.



Operators Manual Section 8.3

- 1. Disinfect as required.
- 2. (A) Remove the Filter Housing and
 - (B) Discard the Filter Element
 - (C) Replace the Filter Housing, but do not insert a replacement Filter at this time. Connect the Transfer Nozzle to the end of the Transfer Line and make sure the Transfer Nozzle Valve is CLOSED.





Operators Manual Section 8.3

- 3. Press the RINSE Start Button. The fill light indicator will turn on and the tank will fill to the 25-Gallon Sensor.
- 4. Once the water reaches the 25-Gallon Sensor the unit will automatically step to RECIRCULATE Operation. Using your safety glasses inspect spray ball operation.
- 5. The unit will step to DRAIN and FILL again.
- 6. When the water reaches the 25-Gallon Sensor during the second RINSE, add 0.5 gallons (1.9 liters) of bleach (sodium hypochlorite 5% to 10%) to the rinse water in the Tank and allow it to recirculate for the duration of the RINSE CYCLE.
- 7. At the completion of the last RINSE Operation (to which the bleach has been added), initiate two (2) complete RINSE CYCLES.

When completed, Check two (2) areas for residual bleach.

(SEE NEXT PAGE)



Operators Manual Section 8.3

Checking from (1) Transfer Hose

- □ Press Dissolution START button. When water has reached the 25 Gallon Sensor, use **STEP MODE** to skip to TRANSFER Operation.
- □ Press the Dissolution START button. Open the Transfer Valve on the Dry Acid Dissolution Unit. Then, slowly open the Valve on the Transfer Nozzle. Allow water to flow to the drain for 30 seconds and then collect a sample to test for residual bleach. Close Transfer Nozzle Valve..

Checking from (2) Drain Hose

- □ Using **STEP MODE** skip to **CYCLE COMPLETE** Operation and press Dissolution **START** Button (132 Gal).
- Using STEP MODE skip to DRAIN Operation and Press Rinse START Button (99 Gal).
- ☐ The Drain Valve will Open. Allow water to drain for 15 seconds and collect a sample from the Drain Hose to test for residual bleach.




Operators Manual Section 8.3

- 8. If residual bleach levels are higher than ANSI/AAMI Standard limit of <0.1 ppm (RD61: 2006), in any of the two areas initiate another complete RINSE Operation. After the RINSE CYCLE is complete, start from section 8.3: number 5, to check for residual bleach. Continue the RINSE CYCLE and test procedure until bleach residuals are within ANSI/AAMI Standard limit of <0.1 ppm (RD61:2006) in both places.
- 9. Once you have attained an acceptable bleach residual reading, connect the Transfer Line back onto the Dry Acid Dissolution Unit Transfer Line Holder.



Operators Manual Section 8.3

- 10. Turn the unit's power OFF and ensure that the main transfer ball valve is CLOSED.
- 11. Remove Filter Housing and drain all residual water from the Housing.
- 12. Install new Filter and tighten Filter Housing into place (See Adjacent Figure).
- 13. Unhook Transfer Nozzle and place Transfer Line onto the Transfer Line Holder.
- 14. Immediately after RINSE CYCLE, make a batch of Dry Acid Product. Leaving the Dry Acid Dissolution Unit with only treated water or wetted with only treated water leaves the Unit susceptible to bacterial growth.





Operators Manual Section 8.3

Sodium Hypochlorite (bleach) Disinfection Step Summary

	Cycle	Operation	Operation #	Comments/Instruction
	Rinse	Fill	1	Press start on the Rinse Cycle
0	Rinse	Recirculate	1	Check for Spray Ball Operation
#	Rinse	Drain	1	
Cycle	Rinse	Fill	2	Press pause and add bleach (0.5 Gal or 1.9 L) then press start to resume.
\leq	Rinse	Recirculate	2	
Ó	Rinse	Drain	2	
	Rinse	Cycle Complete	2	
1	Rinse	Fill	1	Press Start on Rinse Cycle (2 nd Rinse Cycle)
	Rinse	Recirculate	1	
#	Rinse	Drain	1	
Cycle	Rinse	Fill	2	
\Box	Rinse	Recirculate	2	
\sim	Rinse	Drain	2	
<u> </u>	Rinse	Cycle Complete	2	
2	Rinse	Fill	1	Press Start on Rinse Cycle (2 nd Rinse Cycle)
	Rinse	Recirculate	1	
#	Rinse	Drain	1	
<u>_</u>	Rinse	Fill	2	
$\overline{\mathcal{Q}}$	Rinse	Recirculate	2	
Cycle	Rinse	Drain	2	
	Rinse	Cycle Complete	2	
	Dissolution	Fill	Step Mode	
	Dissolution	Transfer	Step Mode	Test from Transfer line**
	Dissolution	Cycle Complete	Step Mode	Test from Drain **

** If residual bleach levels of samples are higher than 0.1 ppm. Continue rinse cycle and test until bleach residuals are less than 0.1 ppm.

After Disinfection is completed,
remember to make a batch of Dry
Acid product.



Transfer Filter and Base Stand Pipe Filter Replacement

Operators Manual: Section 8.4



FILTER REPLACEMENT

Operators Manual Section 8.4

TRANSFER FILTER REPLACEMENT

WHEN?

- The Filter shall be changed AFTER mixing 6 BATCHES.
- When the Dry Acid Dissolution Unit requires DISINFECTION.

NOTE

THE FILTER USED MUST BE COMPATIBLE WITH FRESENIUS MEDICAL CARE DRY ACID PRODUCT AND RATED AT 1 MICRON. FRESENIUS MEDICAL CARE. P/N G84-202-12 IS A POLYPROPYLENE FIBER WOUND ON A POLYPROPYLENE MESH CORE AND MEETS THESE REQUIREMENTS. CELLULOSE FILTERS ARE NOT COMPATIBLE WITH THE FRESENIUS MEDICAL CARE DRY ACID PRODUCT AND WILL BREAK DOWN, CLOGGING AFTER ONLY A FEW BATCHES.




FILTER REPLACEMENT

Operators Manual Section 8.4

TRANSFER FILTER REPLACEMENT

HOW?

- Ensure Mixing Tank is empty
- Power to the Dry Acid Dissolution Unit has been turned off.
- MAIN TRANSFER BALL VALVE is closed.
- Follow the Figure A, B, C, and D for removal and replacement of the filter.





SENSOR REPLACEMENT

Operators Manual: Section 8.5



SENSOR REPLACEMENT

Operators Manual Section 8.5

Sensor Replacement

NOTE: If the final fill sensor needs adjustment or replacement, then a qualified technical personnel shall complete this and the following tasks. Sample must be drawn from the first batch of concentrate Mixed. This sample must be analyzed for correct solution mix before the concentrate can be used. Before replacing the Final Fill Sensor or Relocating the Dry Acid Dissolution Unit contact Fresenius Medical Care Technical Service at 1 (800) 227-2572 and request two (2) empty sample bottles (P/N G83-535-02). See Appendix D of the Operators Manual for further instructions.

(1) Request 2 Empty sample bottles (P/N G83-535-02).



(2) Follow Appendix D

		APPENDIX D: FIRST BATCH VERIFICATION INSTRUCTIONS
1.	Plug in 1	e Final Fill Sensor is replaced or unit relocated and the empty sample bottles are available. the power cord. Turn treated water on. Turn the power switch ON (Red switch on the right he Granufio Dissolution Unit II). Verify the power light activates.
2.	As per ti	ne Operators Manual section 6.0, secure the tank lid and initiate the RINSE CYCLE.
3.	Start Bu	he Operators Manual section 7.0, run a DISSOLUTION CYCLE by pressing the DISSOLUTION tton. When the Add Granules Light flashes check to make sure water has reached the Mid- nsor, and then add the six (6) boxes of GranuFlo® or Citrasate® DRY product.
	①	NOTE: MAKE SURE TO FILL THE PRODUCTION BATCH RECORD FORM ATTACHED TO THIS MANUAL WITH THE PRODUCT INFORMATION.
	\triangle	WARNING: DO NOT USE ANY BOX OF GRANUFLO® OR CITRASATE® DRY CONCENTRATE THAT HAS BEEN OPENED OR TAMPERED WITH. IT IS IMPORTANT THAT THE ENTIRE CONTENTS OF EACH BOX ARE EMPITED INTO THE DISSOLUTION TANK.
	Λ	WARRING! THE USE OF EYE PROTECTION AND GLOVES IS RECOMMENDED WHEN HANDLING DRY ACID PRODUCT. IF CONTACT WITH EYES, RINSE IMMEDIATELY FOR 15 MINUTES. IF CONTACT WITH SKIN, R. LUSH WITH PLEITY OF SOAP AND WATER. SEE MATERIAL SAFETY DATA SHEETS (MSDS) FOR THE DRY ACID PRODUCT BEING USED FOR PURTHER PRESONAL PROTECTIVE EQUIPMENT [POD) OR REMERGINEY REQUIREMENTS, INSTRUCTIONS.
4.		e Operators Manual section 7.1, perform the Specific Gravity Test and record the results on uction batch record form.
5.	analyzed	s sample of the final product(s) using the (2) sample bottles. The product sample(s) will be in accordance to manufacturer's product specifications. If only one product code is being illect two (2) samples of that product to be analyzed.
6.	Place th	e sample bottles and the following completed forms into a shipping box:
	•	Copy of the Production Record Form (Operator's Manual – Form 1 or Form 2)
	•	Batch Analysis Form (Operator's Manual pg 44)
7.	responsi	and affix the pre-printed mailing label from the Batch Analysis Form to the box. It is the bility of the RES to ship the samples to the "Ship To" address on the Batch Analysis Form. Fresenius Medical Care Laboratory at (972) 929-7291 for results.
8.	a qualifi DRY pro	The composition of the first batch of GranuFilo® / Citrasate® DRY product must be tested by ed testing laboratory to ensure that the resulting product meets the GranuFilo® / Citrasate® duct specification. If final solution did not meet final batch criteria for use, the batch of rate must be discarded (See Section 13: Concentrate solution Disposal, Page 37).
	CONCENT	are most be discarded (SEE SECTION 23, CONSENTRATE SOLUTION DISPOSAL, PAGE 37).

(3) Ship Samples Using BATCH ANALYSIS FORM

FRESENIUS MEDICAL CARE
BATCH ANALYSIS FORM Customer Information Form
NOTE: TO PREVENT BATCH ANALYSIS DELAYS, THIS FORM MUST BE COMPLETE AND ACCURATE.
Clinic Name:
Contact Name:Contact Phone Number:
Clinic Fax Number: OR Email: @
Dissolution Tank Serial Number:
Product Catalog Number:
Lot Number:
Date Sample Taken: Sample By:
Fresenius Medical Care Irving Manufacturing 5201 Regent Blvd., Suite 100 Irving, TX 75063 Attention: Laboratory Ship: Overnight



Disposal of Concentrate Solution

Operators Manual: Section 13



Residual Solutions "Bucket Disposal" Section 13.1



"Tank Disposal"
Section 13.2



13.1: RESIDUAL SOLUTION BUCKET DISPOSAL

There are three occasions that require the need to discard residual solution. Those occasions are described as follows:

- **Filter Housing Residual:** Whenever you have to change out a Filter, then the solution in the Filter Housing shall be poured into the *Residual Solution Bucket*.
- Hydrometer Cylinder Residual: Once you are complete with the Specific Gravity Test, the solution in the Hydrometer Cylinder shall be poured into the Residual Solution Bucket.
- Initial 3.5 Gallon Transferred Solution: This is the 3.5 gallons of solution transferred into the *Residual Solution Bucket* in the beginning part of the Specific Gravity Test, section 7.1.



13.2: TANK SOLUTION DISPOSAL

There are five occasions in which you will need to discard the solution in the Dry Acid Dissolution Unit. Those occasions are described as follows:

- **Expired Solution:** If solution remains in the Dry Acid Dissolution Unit Tank for more than 14 days.
- Foreign Object: Any foreign object falling into the Dry Acid Dissolution Unit Tank after the dry acid product has been added.
- Incorrect Catalog Used: Any batch preparation with incorrect dry acid product catalog (s) number (s).
- Specific Gravity Out of Range: Any batch prepared that has been found out of range.
- **Dry Acid Unit Needing Service:** Any Dry Acid Dissolution Unit that contains solution that becomes inoperable.



continued

CONCENTRATE SOLUTION DISPOSAL PROCEDURE:

MIXED DRY ACID CONCENTRATE SOLUTION HAS A PH BELOW 6. THEREFORE, IF YOU MUST DISPOSE OF ACID CONCENTRATE SOLUTION BECAUSE OF ONE OF THE ABOVE SITUATIONS, YOU SHOULD FOLLOW ANY APPLICABLE DISPOSAL REQUIREMENTS OF YOUR LOCAL, STATE, AND/OR FEDERAL AUTHORITIES. SEE CLINIC MANAGER FOR MORE INFORMATION.



WARNING! Some Chemicals used to neutralize acid concentrate solutions may cause splattering and/or generate dangerous levels of gases when combined. For example, potentially dangerous levels of carbon dioxide may be released when sodium bicarbonate is used as the neutralizing agent. Please carefully consider these issues, including proper ventilation, if neutralization is required.



continued

Immediately after a tank disposal, step unit to CYCLE COMPLETE, press STEP MODE button to turn OFF, then ensure that the water inlet valve is open and perform two (2) RINSE CYCLES by pressing the RINSE CYCLE start button and by pressing it again when the first RINSE CYCLE is completed.

After completing two (2) full RINSE CYCLES visually inspect if there is any powder deposits remaining in the tank. These can be removed by running additional RINSE CYCLES in the unit.


Immediately after rinsing the tank clean, make a batch of Dry Acid Product. Leaving the Dry Acid Dissolution Unit with only treated water or wetted with only treated water leaves the Unit susceptible to bacterial growth.

CAUTION! DO NOT ALLOW THE UNIT TO REMAIN FULL OF WATER WITHOUT THE ADDITION OF FRESENIUS MEDICAL CARE DRY ACID PRODUCT. BACTERIAL GROWTH MAY OCCUR.




Manual Control Operations












Pause State

PAUSE STATE

If at any time the operator needs to PAUSE a timed operation during the cycle, the PAUSE button may be pressed.

- This will cause the Indicating light for the current step of the operation to flash.
- In the PAUSE state Pressing the PAUSE button or placing the control into the STEP MODE will disable the pump, agitators, drains, fill valves, etc.
- To continue the cycle, press the START button and the timers will resume.






Step Mode


When to Use it?

The STEP MODE function is intended to be used during the Disinfection Operation or when it is necessary to discard an incorrectly mixed batch of solution.

How To Use it?

- Press the STEP MODE ON/OFF button
- The step Mode Indicator light will illuminate.
- System will enter STEP MODE and all Operations will be suspended.
- Press STEP button to skip to desired Operation.
 (operation indicator light will illuminate)
- Press the STEP MODE ON/OFF button and the Operation is continued.







Main Power Disconnect

MAIN POWER DISCONNECT

The MAIN POWER Switch is provided to allow the operator to completely shut down the power to the Dry Acid Dissolution Unit.

The Main Power Disconnect Switch should be switched to the OFF position when the Dry Acid Dissolution Unit is not in use or in case of an emergency.

Remove POWER PLUG from wall receptacle to disconnect power. A 'LOCKOUT' device may be used to prevent unauthorized start up.



Appendix and Forms



OPERATORS MANUAL APPENDIX A

SPECIFIC GRAVITY RANGES

APPENDIX A: GRANUFLO® DRY ACID SPECIFIC GRAVITIES TABLE

Specific Gravity Ranges

1:44 PROPORTIONING

		17°C to 21°C		22°C to 26°C		27°C to 31°C	
Catalog Number		Low	High	Low	High	Low	High
1K	OFD1251-3B	1.191	1.203	1.188	1.200	1.186	1.198
	0FD2201-3B	1.191	1.203	1.189	1.201	1.187	1.199
2K	OFD2231-3B	1.192	1.204	1.190	1.202	1.187	1.199
7	0FD2251-3B	1.192	1.204	1.190	1.202	1.188	1.200
	0FD2301-3B	1.193	1.205	1.191	1.203	1.189	1.201
	0FD3201-3B	1.193	1.205	1.191	1.203	1.188	1.200
3K	0FD3231-3B	1.194	1.206	1.191	1.203	1.189	1.201
9	0FD3251-3B	1.194	1.206	1.192	1.204	1.189	1.201
	OFD3301-3B	1.195	1.207	1.193	1.205	1.190	1.202

Minimum Thermometer Specifications: Temperature range 25° C +/-5°C (68° to 86° F) and accuracy +/- 1°C (3.6 °F)

For reference only.

See Appendix A in Operator's Manual



OPERATORS MANUAL APPENDIX B

SPECIFIC GRAVITY RANGES

APPENDIX B: CITRASATE® DRY SPECIFIC GRAVITIES TABLE

Specific Gravity Ranges

1:44 PROPORTIONING


	16.5°C to 21.4°C (61.7°F to 70.6°F)		21.5°C to (70.7°F to		26.5°C to 31.4°C (79.7°F to 88.5°F)		
Catalog Number		Low	High	Low	High	Low	High
	0FD1231-DA	1.192	1.204	1.189	1.201	1.187	1.199
1K	0FD1251-DA	1.192	1.204	1.190	1.202	1.188	1.200
	0FD2231-DA	1.194	1.206	1.191	1.203	1.189	1.201
X	0FD2251-DA	1.194	1.206	1.192	1.204	1.189	1.201
7	0FD2301-DA	1.195	1.207	1.193	1.205	1.190	1.202
	0FD3231-DA	1.195	1.207	1.193	1.205	1.191	1.203
3	0FD3251-DA	1.196	1.208	1.194	1.206	1.191	1.203
_ c	0FD3301-DA	1.197	1.209	1.194	1.206	1.192	1.204


Minimum Thermometer Specifications: Temperature Range 25° C +/-5°C (68° to 86° F) and accuracy +/1°C (3.6°F)

For reference only.
See Appendix B in Operator's Manual



OPERATORS MANUAL APPENDIX C







OPERATORS MANUAL APPENDIX D

For reference only. See Appendix D in Operator's Manual

APPENDIX D: FIRST BATCH VERIFICATION INSTRUCTIONS

- Once the Final Fill Sensor is replaced or unit relocated and the empty sample bottles are available. Plug in the power cord. Turn treated water on. Turn the power switch ON (Red switch on the right side of the Dry Acid Dissolution Unit). Verify the power light activates.
- As per the Operators Manual P/N 460018 section 6.0, secure the tank lid and initiate the RINSE CYCLE.
- As per the Operators Manual P/N 460018 section 7.0, run a DISSOLUTION CYCLE by pressing the DISSOLUTION Start Button. When the ADD GRANULES light flashes check to make sure water has reached the Mid-Level Sensor, and then add the eight (8) boxes of GranuFlo® or Citrasate® DRY product.



NOTE: MAKE SURE TO FILL THE PRODUCTION BATCH RECORD FORM ATTACHED TO THIS MANUAL WITH THE PRODUCT INFORMATION.



WARNING: DO NOT USE ANY BOX OF GRANUFLO® OR CITRASATE® DRY CONCENTRATE THAT HAS BEEN OPENED OR TAMPERED WITH. IT IS IMPORTANT THAT THE ENTIRE CONTENTS OF EACH BOX ARE EMPTIED INTO THE DISSOLUTION TANK.



WARNING: USE PROTECTIVE EYE WEAR AND LATEX (OR EQUIVALENT) GLOVESWHEN HANDLING GRANUFLO®.

- As per the Operators Manual P/N 460018 section 7.1, perform the Specific Gravity Test and record the results on the production batch record form.
- Collect a sample of the final product(s) using the (2) sample bottles (P/N G83-535-02). The product sample(s) will be analyzed in accordance to manufacturer's product specifications. If only one product code is being used, collect two (2) samples of that product to be analyzed.
- 6. Place the sample bottles and the following completed forms into a shipping box:
 - Copy of Batch Production Record Form (Operator's Manual Form 1 or Form 2)
 - Batch Analysis Form (Operator's Manual pg 43)
- Cutout and affix the pre-printed mailing label from the Batch Analysis Form to the box. It is the
 responsibility of the RES to ship the samples to the "Ship To" address on the Batch Analysis Form.
 Contact Fresenius Medical Care Laboratory: at (972)-929-7291 for results.
- 8. NOTICE: The composition of the first batch of GranuFlo® / Citrasate Dry product must be tested by a qualified testing laboratory to ensure that the resulting product meets the GranuFlo® / Citrasate Dry product specifications. If final solution did not meet final batch criteria for use, the batch of concentrate must be discarded (SEE SECTION 13: CONCENTRATE SOLUTION DISPOSAL PROCEDURES, PAGE 37).



OPERATORS MANUAL

FORM 1: GranuFlo® Batch Production Record

FORM 1: GRANUFLO® BATCH PRODUCTION RECORD

Dialysis Unit Name & L	ocation #:			Ory Acid Dissolution (Jnit Serial #:	
DRY ACID PRODUCT CAS	EINFORMATION					
OPERATOR (print name):		DATE:		1	TIME:	Batch #:
CASE 1: Dry Acid Catalog # (label on box)			BOX Lot #			Potassium # (1K, 2K, 3K)
CASE 2: Dry Acid Catalog #		E	BOX Lot #			Potassium #
CASE 3: Dry Acid Catalog #		E	BOX Lot #			Potassium #
CASE 4: Dry Acid Catalog #			BOX Lot #		Potassium #	
CASE 5: Dry Acid Catalog #			BOX Lot #		Potassium #	
CASE 6: Dry Acid Catalog #		E	BOX Lot #		Potassium #	
CASE 7: Dry Acid Catalog #		E	BOX Lot #			Potassium #
CASE 8: Dry Acid Catalog #		E	BOX Lot #			Potassium #
1. After Final Fill Level ha	s been reached, CLOSE water valve to t	the Unit. Once				
Measured Temp	Print Catalog # -Specific Gravity Value	for the Measured	SPECIFIC GR		Measured Specific Gravity Value:	Check one
TEMP:	LOW:	HIG		represent A	measured opening Gravity Value.	□-Pass
						☐-Fail (void section)
OPERATOR SIGNATURE:	<u>'</u>			VERIFIER SIGNATU	RE:	, , ,

For reference only.
See form in Operator's Manual



OPERATORS MANUAL

FORM 2: Citrasate® DRY Batch Production Record

DIALYSIS UNIT NAME &	LOCATION#:		Dry Acid Dissolution U	Init Serial #:		
DRY ACID PRODUCT CA	ASE INFORMATION					
OPERATOR (print name):		DATE:	1	TIME:		Batch #:
CASE 1: Dry Acid Catalog #	(label on box)	BOX Lot #			4 [™] bag used ✓ box	Potassium # (1K, 2K, 3K)
					Yes 🗆	
CASE 2: Dry Acid Catalog #		BOX Lot #			4 [™] bag used ✓ box	Potassium #
					Yes 🗆	
CASE 3: Dry Acid Catalog #		BOX Lot #			4™ bag used ✓ box	Potassium #
					Yes 🗆	
CASE 4: Dry Acid Catalog #		BOX Lot #			4 [™] bag used ✓ box	Potassium #
					Yes 🗆	
CASE 5: Dry Acid Catalog #		BOX Lot #			4™ bag used ✓ box	Potassium #
Name and Address of the Control of t					Yes 🗆	
CASE 6: Dry Acid Catalog #		BOX Lot #			4 [™] bag used ✓ box	Potassium #
					Yes 🗆	
CASE 7: Dry Acid Catalog #		BOX Lot #			4 th bag used ✓ box	Potassium #
					Yes 🗆	
CASE 8: Dry Acid Catalog #		BOX Lot #			4 [™] bag used √ box	Potassium #
					Yes 🗆	
			SPECIFIC GRAVITY			
Measured Temp		lic Gravity Value for the Measure		Measured Sp	pecific Gravity Value:	Check one
TEMP:	LOW:	HIGH:	1			□-Pass
						☐-Fail (void section)
1. After Final Fill Level	has been reached, CLOSE water	r valve to the Unit. Once th	nis is done check the box			
OPERATOR SIGNATURE:			VERIFIER SIGNATUR	RE:		

For reference only.
See form in Operator's Manual




BASE STAND PIPE FILTER

(Only on the 99 Gal, Operators Manual Sec 8.4.2)

BASE STAND PIPE FILTER REMOVAL

- Using Base Filter Rod, P/N 260025 (A)
- Reach into the tank and connect the small end of the shaft into the Screw Nut of Base Filter (B).
- Turn the Screw Nut Counter Clockwise until the Base Stand Pipe Filter is no longer attached to the tank.
- Use the Base Filter Rod to help lift the Base Stand Pipe Filter out of Tank.



After Removing the Base Stand Pipe Filter.

- Clean with purified water
- Re-install



Copyright and Trademark Statements

© 2013, 2016, 2017, 2023 Fresenius Medical Care. All Rights Reserved. Fresenius Medical Care, the triangle logo, GranuFlo and Citrasate are trademarks of Fresenius Medical Care Holdings, Inc. or its affiliated companies. All other trademarks are the property of their respective owners.

